

INSIDE FOR TRAINING



INSIDE from Faro is the first Endoscope made specifically for Dentistry. Its fundamental characteristic is its ability to manage images intelligently.

For this reason, the **INSIDE is also an excellent tool for training establishments and laboratories.**


Current solutions involve either ceiling cameras where the trainers' head often gets in the way or hand held intra-oral cameras that have poor resolution and require an assistant to hold it steady to avoid it being knocked out of place.

The **INSIDE** system however, can be easily attached to existing training equipments or be free standing / ceiling mounted as required for convenience. *The cantilevered arms ensure the INSIDE camera head remains stable whilst the trainer has both hands free to work.* The images can be rotated and shots captured with ease from the control panel.

Images from the **INSIDE** can then be sent out to any number of other screens and monitors, (e.g. student training positions and / or to a central projector / LCD screen) allowing the trainer to show procedures in superb resolution, highly magnified and without the annoying issues mentioned earlier.

Furthermore, the system can be set up to allow the trainer to record training procedures to digital recorders (for future playback or use with on-line libraries) as well as show pre-saved images from a computer. Touch screen monitors would allow the trainer to annotate over images in real time thus further enhancing the learning experience for students.

Finally, not only can the system be used in one central location, but the trainer could actually train or offer specialist advice between locations using simple video and data conferencing techniques.



To see INSIDE in action visit
www.phoenix-endoscope.co.uk
Or call us on
01939 200467

